

enovis™

*DONJOY*®

MaxTrax® Xcel ROM

**SIMPLE, RELIABLE AND PRACTICAL.**



# MaxTrax<sup>®</sup> Xcel ROM

Simple, reliable and practical.

When you need a **cost-effective range-of-motion (ROM) walker boot**, look no further than MaxTrax<sup>®</sup> Xcel ROM. Intended to aid recovery from a variety of ankle and foot indications, its **tried-and-tested hinge allows safe ankle immobilization as well as ROM control** between 45° plantar-flexion and 30° dorsi-flexion. Together with its shock-absorbing rocker sole and padded liners, this is how MaxTrax Xcel ROM provides **comfortable and supportive ROM control at great value.**

## FEATURES AND BENEFITS:

- Adjustable ROM hinge: **Allows ROM settings in 7.5° increments between 45° plantar-flexion and 30° dorsi-flexion** with locking option
- Unique pre-shaped double upright design: **Modifies pressure** on the malleoli and lower leg and **ensures a better fit**
- Edgeless rocker sole: Helps promote a **natural, stable gait**
- Shock-absorbing sole: **Reduces the impact of heel-strike** during ambulation
- Padded removable liners: For cool **comfort** and **easy cleaning and replacement**
- Adjustable hook-and-loop straps: Help **reduce pressure** on the injury site
- XS-XL sizes, tall and ankle versions: Accommodate a **range of leg sizes and heights**



## INTENDED USE/INDICATIONS:

- Post-operative or post-traumatic immobilization of the ankle or foot
- Severe ankle and/or foot sprain
- Conservative treatment of ankle and/or foot fractures
- Severe ankle tendinopathy
- Conservative treatment of Achilles tendon rupture

## ORDERING INFORMATION

PART NO.		EU SHOE SIZE	SIZE
ANKLE	TALL		
82-102-XS	82-104-XS	36	XS
82-102-S	82-104-S	37 - 40	S
82-102-M	82-104-M	41 - 44	M
82-102-L	82-104-L	45 - 47	L
82-102-XL	82-104-XL	48+	XL

**enovis**<sup>™</sup>

Creating Better Together<sup>™</sup>

**enovis-medtech.eu**

1a Guildford Business Park, Guildford, Surrey, GU2 8XG, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice.

The contents of this document do not constitute medical, legal, or any other type of professional advice.

Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-5618245833-EN - Rev A

Copyright © 2023 DJO, LLC, a subsidiary of Enovis Corporation